

## Product Information

Date: 3/1//2010

by:KRD

### Medium Capacity Wedge Grips

<p><b>TH243-20 Wedge Grip</b></p> <p>Opens to 13 mm (0.512")          Face dimensions 35 mm x 50 mm (1.38" x 1.97")          Pyramid jaws face standard, serrated available          Lever on gear rack assist for closing and opening          Rated to 20 kN (4600 lbf)          Available in pairs only</p>	
<p><b>GF-9 Wedge Grip</b></p> <p>Opens to 6mm (0.25")          Face dimensions 25mm x 35mm (0.98" x 1.38")          Serrated jaw face only          Direct acting lever assist for opening and closing          Rated to 2.5 kN (565 lbf)          Available in pairs only</p>	
<p><b>432-398 Wedge Grip</b></p> <p>Opens to 8mm ( 0.314")          Face dimensions 25 mm x 30 mm (0.98 x 1.18")          Serrated jaw face only          Direct acting lever assist for opening and closing          Rated to 5 kN (1100 lbf)          Available individually</p>	
<p><b>ML-1700 Miniature Tensile Grip</b></p> <p>Opens to 5.56 mm (0.22")          Face dimensions 6mm x 25mm (0.236" x 0.98")          Serrated jaw face only          Direct acting lever assist for opening and closing          Rated to 2 kN (450 lbf)          Available individually</p>	
<p><b>TG15 Wedge Grip</b></p> <p>Opens to 10 mm (0.39")          Face dimensions 25mm x 50mm (0.98" x 1.97")          Cross cut jaw face standard, V notch and Straight cut avail          Direct acting lever assist for opening and closing          Rated to 5 kN (1100 lbf)          Available in pairs only</p>	